



Part C. DNase Digestion and clean-up [microRNA Profiling]

Adjust Dnase digestion volume to 50 ul if you have less then 15 ug of LMW RNA.

- 1 Dilute 30 ug of LMW RNA in a volume of 84 ul of Qiagen RNase-free water.
- 2 Add 10ul of Invitrogen 10X Dnase Digestion Buffer
- 3 Add 2ul of Rnase Block (Rnase Inhibitor) -mix by gentle vortexing
- 4 Add 4ul of Ambion Rnase-free Dnase I
- 5 Mix by inversion or flicking and spin briefly (do not vortex)
- 6 Incubate at 37C for 30 min
- 7 Add 4.0ul of 0.5 M EDTA pH 8 , vortex and spin down
- 8 Adjust volume of Dnased-RNA to 100 ul with Qiagen water if necessary
- 9 Combine DNased digest RNA with 350 ul of buffer RLT and mix by gentle vortexing (don't centrifuge)
- 10 Combine RNA/RLT mixture with 1500 ul of 100% ethanol
- 11 Mix by gentle vortexing (don't centrifuge)
- 12 Apply the sample to an Rneasy MinElute spin column in a 2ml collection tube.
- 13 Close the tube gently and centrifuge for 1 min at >10,000 g.
(Repeat till all the sample is loaded on to the column.)
- 14 Transfer the spin column into a new 2ml collection tube.
- 15 Pipet 500 ul buffer RPE onto the spin column
(Apply around rim of column to make sure the area is washed)
- 16 Close the tube gently, and centrifuge for 1 min at >10,000 g to wash the column.
- 17 Discard the flow-through and reuse collection tube
- 18 Pipette 500 ul of 80% ethanol to the Rneasy minelute spin column.
- 19 Close the tube gently, and centrifuge for 1 min at >10,000 g
- 20 Discard flow-through and collection tube.
- 21 Transfer the Rneasy minelute spin column into a new 2 ml collection tube.
- 22 Close the cap of the spin column, and centrifuge in a microcentrifuge at 16,000 g or max speed for 3 min.
- 23 Discard flow-through and collection tube.
- 24 Repeat steps 21-23 one more times (important to remove any trace of ethanol)
- 25 To elute, transfer the spin column to a new 2 ml capless collection tube.
- 26 Pipette 14 ul Rnase-free water directly onto the center of the silica-gel membrane.
- 27 Close the tube gently, incubate 1 minute on bench and centrifuge for 1 min at max speed to elute
- 28 Measure the sample concentration using Ribogreen based fluorescence assay. (optional)
- 29 Freeze at -70C or ship to ORB on dry ice.